
# THE FRANCK-HERTZ EXPERIMENT KITTY HARRIS AND JOSH ELSARBOUKH

### ELECTRON ENERGY AND SCATTERING CROSS-SECTION




Cross-section of inelastic collisions vs. electron energy. Uncertainty of 30% due to measurement imprecisions. Haken & Wolf, p 698

The Franck-Hertz experiment, presented in 1914 by James Franck and Gustav Hertz, involves accelerating electrons through mercury (also done with neon) vapor via a potential grid, towards a beam current collecting anode. As the accelerated electrons gain sufficient energy for excitation, the cross-section of inelastic collisions increases. This can be seen in the diagram to the left, which shows roughly the cross section of inelastic electron-mercury collisions as a function of electron energy. The points of highest cross-section correspond approximately to the three lowest energy states at 4.67, 4.89, and 5.46 volts, respectively.

The cross-section can best be understood as an effective 'area' inside of which a scattering event can occur. It is not only related to accelerating electron energy, but also the mean free path of the electrons, as well as temperature.

$$\sigma = \frac{k_B T}{P \lambda}$$


(Where P is the pressure, which is a function of the temperature T, and  $\lambda$  is the mean free path, which is dependent on electron energy.)



Lowest energy states of mercury. Haken & Wolf, p 697

# CURRENT MINIMA AT ENERGY INTERVALS

 To reach the collection plate, electrons must overcome a potential difference.



• Only electrons with sufficient energy can overcome this difference:

$$\frac{m_e}{2}v^2 = KE \ge U = q\Delta V$$

So electrons which have lost too much energy in collisions can't make it across the 1.5[V] potential.

- Current reflects this: I = neAv increases as the number of electrons (that reach the collection plate, where the current is being read) increases.
- We see a smooth curve because velocity (and therefore energy) is a probabilistic
   distribution, not the same for each electron, and because not all electrons will collide.

# **SETUP**


We will be focusing on the Keithley 6487 and the portion inside the quartz tube.

Diagrams uploaded by Dr. Tagg to

https://sites.google.com/site/experimentalphysicsdecathlon/home/
04-fundamental-quantumbehavior/franck-hertzexperiment.

Diagram to the right is titled "FRANCK HERTZ.bmp".

Diagram in the next slide is pulled from the safety instructions in "device.pdf".



# INSIDE THE QUARTZ TUBE

#### Quartz Tube

- Electromagnetically Non-Conductive
- Transparent
- Common: Makes it relatively cheap compared to alternatives
- Average thermal conductivity prevents it from having a major contribution to temperature.
- Vacuum prevents vapor from escaping, prevents other gases from entering, and keeps pressure low.

#### Mercury Vapor

- Boiling Point: 629.88[K] much greater than the box.
  - However, pressure is very low  $(O\sim 10^{-3} [atm])$ , so vapor is easy to maintain.
- Excited states are distinct.


#### Tungsten Fillament (Cathode)

- Electrons are boiled off high resistance > excess heat > high electron energy
   Grid (Anode)
- This is where we sweep potential relative to the cathode.

#### Collection Plate

- Swept to stay at a 1.5[V] difference from the anode, reducing KE of incoming electrons.
- This is where current is measured.





# KEITHLEY 6487



- Ammeter Resolution: 10[fA]
- Sweeping Voltage:
  - 200 microvolts to 500 volts
  - 200 steps per second
  - 0.2 millivolt resolution
- Burden Voltage: <200 microvolts

- Dual picoammeter and variable voltage source
- Ideal for low-, dark-, and beam-current measurements
- Fine resolution, low noise
- I/O trigger modes for automated production

Range	5½ Digit Default Resolution	Accuracy (1 Year) <sup>1</sup> ±(% rdg. + offset) 18°–28°C, 0–70% RH	Typical RMS Noise <sup>2</sup>	Typical Analog Rise Time (10% to 90%)³ Damping⁴ Off On	
2 nA	10 fA	0.3 % + 400  fA	20 fA	4 ms	80 ms
20 nA	100 fA	0.2 % + 1 pA	20 fA	4 ms	80 ms
200 nA	1 pA	0.15% + 10  pA	1 pA	$300 \mu \mathrm{s}$	1 ms
2 μΑ	10 pA	0.15% + 100 pA	1 pA	$300 \mu s$	1 ms
20 μΑ	100 pA	0.1 % + 1  nA	100 pA	$110 \mu s$	$110 \mu s$
200 μΑ	1 nA	0.1 % + 10  nA	100 pA	110 μs	110 μs
2 mA	10 nA	0.1 % + 100  nA	10 nA	$110 \mu s$	$110 \mu s$
20 mA	100 nA	$0.1 \% + 1 \mu A$	10 nA	110 μs	110 µs

# TEMPERATURE-DEPENDENCE OF ΔE

Temperature [K]	433	447	462	476
Energy (E <sub>a</sub> ) [eV]	4.85	5.03	4.99	4.91

Uncertainties are  $\pm 0.1$ [eV]

Since our carriers are electrons, voltage [V] and energy [eV] measurements are identical.


Minima were determined by taking curve fits to slices of data.

Then the distances between these minima were plotted.

$$\Delta E = \left[1 + \frac{\lambda}{L}(2n - 1)\right] E_a$$

$$\sigma = \frac{k_B T}{p \lambda}$$

$$E_a = \Delta E_{(n = \frac{1}{2})}$$




# FINDING THE CROSS-SECTION FOR COLLISIONS

 Cross-section is sensitive to movement of particles, therefore varying with temperature:

$$\sigma = \frac{k_B T}{8.7 * 10^{9 - \frac{3110[K]}{T}} [Pa] \lambda}$$

• The mean free path is taken from the slope of the function:  $\frac{d\Delta E}{dn} = \frac{2\lambda}{L} E_a$ 

Temperature [K]:	433	447	462	476
λ [10 <sup>-6</sup> m]:	41.2	1.59	6.41	12.2
Uncertainty:	9.24	8.90	8.96	9.12
$\sigma$ [10 <sup>-18</sup> m $^2$ ]:	2.54	40.5	6.16	2.11
Uncertainty [m²]:	0.569	22.7	8.62	1.58



# **WORKS CITED**

- Abdulagatov, I. M., et al. "Thermal Conductivity of Fused Quartz and Quartz Ceramic at High Temperatures and High Pressures." ScienceDirect, Academic Press, May 2000, www.sciencedirect.com/science/article/pii/S0022369799002681.
- "Alpha: Making the World's Knowledge Computable." Wolfram, 2019, www.wolframalpha.com/.
- Friedman, Hershel. "Minerals.net." Ice: The Mineral Ice Information and Pictures, 2019, www.minerals.net/mineral/quartz.aspx.
- <u>Hackman. "Information on Electronic Quartz Crystals." Hackman's Realm, lateblt.tripod.com/xtals.htm.</u>
- Hann, G.F. "What Really Happens in the Franck-Hertz Experiment with Mercury?" American Journal of Physics 56, 696 (1988). http://home.uni-leipzig.de/~physik/sites/mona/wp-content/uploads/sites/3/2017/08/Franck-Hertz-Experiment.pdf.
- "What Is Important About Thermal Properties?" Amethyst Galleries Mineral Gallery, 2014, www.galleries.com/minerals/property/thermal.htm.
- "5: Thermionic Emission." http://www.physics.csbsju.edu/370/thermionic.pdf.