
DIFFRACTION

Experimental Collaboration From University of Colorado Denver

Kathryn Harris and Idriss Kacou

Associate Professor: Dr. Carlson





- Introduction
- Core Ideas
 - The Diffraction Pattern
 - The wave aspects of light
 - Effects of aperture width
 - Optical geometry and optical wavelength
 - The location of peaks
 - Fresnel (or Near) Field and Fraunhofer (or Far) Field
 - Other Aspects of the Light Source
 - Wavelength
 - The light source may effect a diffraction pattern
 - The aperture and the coherence of the source
 - The light intensity
- Experiment
 - CCDs
 - Laser
- Data Analysis
 - Data fitting to the curves or Newton's method
 - Determining the angles of diffracted light without damaging the CCD

CORE IDEAS

The diffraction of the light source



THE DIFFRACTION PATTERN


- The wave aspects of light
- Effects of aperture width
- Optical geometry and optical wavelength
- The location of peaks

$$\frac{D}{2}\sin\theta = \frac{\lambda}{2}$$

 Fresnel or near field and Fraunhofer or far field


$$F = \frac{R^2}{D\lambda}$$







- Wavelength
 - $y = A \sin(kx \omega t)$
- The light source may effect a diffraction pattern
- The aperture and the coherence of the source
- The light intensity
 - Beam Quality

EQUIPMENT










What is not listed here is our computer setup. The CCD we used was a webcam with the lens removed; this plugged directly into the computer, which we used to operate it and save our data. This data was then translated into a matrix by MatLab, which we used to do our data analysis.

Distances were measured externally to avoid damaging any of the equipment.


EQUIPMENT AND EXPERIMENTAL SETUP

CCDS LASERS






- Measure intensity of light
- Multiple electrodes per pixel measure different colors
- Potential wells attract electrons more electrons at brighter regions
- Quantum Efficiency
- Noise Physical damage, Thermal



- Source which produces a very narrow beam of light
- It is different from a light bulb or a flash light
- Beam Quality
 - $I(r) = I_0 e^{-2(\frac{r}{w})^2}$
 - I= peak intensity (in the center of spot)
 - W= beam radius (1/e^2)



DATA ANALYSIS


CALCULATED PIXEL SIZES USING ANGLES

Slit:/Wavelength:	650[nm] (R)	532[nm] (G)	405[nm] (V)
0.04[mm]	18[μm]±6[μm]	29[μm]±10[μm]	1.6[μm]±0.2[μm]
0.08[mm]	11[μm]±3[μm]	6.4[μm]±2.1[μm]	1.6[μm]±0.1[μm]
0.16[mm]	1.8[μm]±0.6[μm]	2.4[μm]±0.9[μm]	2.1[μm]±0.5[μm]

$$\frac{D}{2}\sin\theta = n\frac{\lambda}{2}$$
 Substitutes into $\frac{ns}{d} = \sin\theta$ Which we can average to get our pixel size.





V 0.04





G 0.04







REFERENCES

- CCDs:
 - McFee, Chris. "An Introduction to CCD Operation." UCL Department of Space & Climate Physics Mullard Space Science Laboratory, www.mssl.ucl.ac.uk/www_detector/ccdgroup/optheory/ccdoperation.html
- Near and Far Fields:
 - Laserist. "Optics Notes." Near Field and Far Field, 1 Jan. 1970,
 laseristblog.blogspot.com/2011/12/near-field-and-far-field.html

